



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Oct-1/2 Monoclonal Antibody

单克隆抗体

产品货号: MA01797

第三版

免疫原:Synthetic Peptide

亚型:IgG

克隆:5G3

描述:The OCT1 transcription factor was among the first identified members of the POU transcription factor family. Members of this family contain the POU domain, a 160-amino acid region necessary for DNA binding to the octameric sequence ATGCAAAT.

配方:

浓度:1mg/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF