



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human A1BG Monoclonal Antibody

单克隆抗体

产品货号: MA01837

第三版

免疫原:Fusion protein of A1BG

亚型:IgG1

克隆:507

描述:A1BG (alpha-1B-glycoprotein) is a plasma protein and is a member of the immunoglobulin superfamily. A1BG contains 474 amino acids and has 5 intrachain disulfide bonds. It has been reported that cysteine-rich secretory protein 3 (CRISP-3) is a specific and high-affinity ligand of A1BG. A1BG-CRISP-3 complex displays a similar function in protecting the circulation from a potentially harmful effect of free CRISP3 .

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF