

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PSAP Monoclonal Antibody

单克隆抗体

产品货号: MA01839

第三版

免疫原:Fusion protein of ACPP

亚型:IgG2a

克隆:220

描述:ACPP is also named as PAP(prostatic acid phosphatase),5'-nucleotidase,thiamine monophosphatas and Belongs to the histidine acid phosphatase family.This enzyme is a homodimeric glycoprotein with a subunit molecular mass of 50 kDa and is an androgen-stimulated secretory product of the human prostate gland.It is also a non-specific phosphatase and has been implicated in the progression of prostate cancer,which inactivates LPA in human seminal plasma. This protein has three glycosylation sites.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF