

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ADIPOQ Monoclonal Antibody

单克隆抗体

产品货号: MA01840

第三版

免疫原:Fusion protein of ADIPOQ

亚型:IgG1

克隆:588

描述:Adiponectin (AdipoQ),an adipocyte-derived hormone,is one of the most abundant adipokines in the blood circulation. Adiponectin modulates a number of metabolic processes,including improving insulin sensitivity and anti-inflammatory activity. The role of AdipoQ in reproduction is not yet fully understood,but the expression of AdipoQ in reproductive tissues has been observed in various animals and humans,including chicken testis,bovine ovary,and human placenta. Adiponectin exerts its effects by activating a range of different signaling molecules via binding to two transmembrane AdipoQ receptors,AdipoR1 and AdipoR2. AdipoR1 is expressed primarily in the skeletal muscle,whereas AdipoR2 is predominantly expressed in the liver. AdipoQ May play a role in cell growth,angiogenesis and tissue remodeling by binding and sequestering various growth factors.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF