
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse AFM Monoclonal Antibody

单克隆抗体

产品货号: MA01841

第三版

免疫原:Fusion protein of AFM

亚型:IgG1

克隆:790

描述:AFM (Afinin) is a member of the albumin superfamily,which comprises albumin,vitamin D-binding protein, α -fetoprotein,and afamin. AFM is present in the plasma/serum,cerebrospinal fluid,and follicular fluid,and functions as a vitamin E-binding protein. Alterations of expression level of AFM have been linked to variety of diseases like cancer.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF