

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human AHCY Monoclonal Antibody

单克隆抗体

产品货号: MA01842

第三版

免疫原:Fusion protein of AHCY

亚型:IgG2b

克隆:190

描述:AHCY,also named as SAHH,belongs to the adenosylhomocysteinase family. AHCY is a competitive inhibitor of S-adenosyl-L-methionine-dependent methyl transferase reactions; therefore AHCY may play a key role in the control of methylations via regulation of the intracellular concentration of adenosylhomocysteine. AHCY catalyzes the reversible hydrolysis of S-adenosylhomocysteine to adenosine and L-homocysteine..

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF