

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human AK2 Monoclonal Antibody

单克隆抗体

产品货号: MA01843

第三版

免疫原:Fusion protein of AK2

亚型:IgG2a

克隆:205

描述:AK2(Adenylate kinase 2,mitochondrial) is also named as ADK2 and belongs to the adenylate kinase family. It plays a unique role in energy metabolism and energy transfer by regulating the ATP/ADP rate between the cytoplasm matrix and the mitochondria and it is increased expression during adipocyte and B cell differentiation. AK2 catalyzes the reversible transfer of a phosphoryl group from ATP to AMP. Defects in AK2 are the cause of reticular dysgenesis (RDYS). It has 6 isoforms produced by alternative splicing.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF