



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat CD166 Monoclonal Antibody

单克隆抗体

产品货号: MA01846

第三版

免疫原:Fusion protein of Albumin

亚型:IgG1

克隆:384

描述:Albumin is the most abundant protein in blood plasma. Alterations of level of serum albumin are linked to variety of diseases. Albumin is expressed exclusively by well-differentiated hepatocytes,thus anti-albumin has been used to mark hepatocytes.
(21388516,23832071)

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF