

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human AMOT Monoclonal Antibody

单克隆抗体

产品货号: MA01853

第三版

免疫原:Synthetic Peptide

亚型:IgG2b

克隆:678

描述:predominantly in endothelial cells of capillaries as well as larger vessels of the placenta where it may mediate the inhibitory effect of angiostatin on tube formation and the migration of endothelial cells toward growth factors during the formation of new blood vessels. Most abundant expression was found in placenta and skeletal muscle. AMOT has two isoforms with MW 130 kDa (p130) and 80 kDa (p80). This antibody recognizes p130 only.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF