

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat AMY2A Monoclonal Antibody

单克隆抗体

产品货号: MA01854

第三版

免疫原:Fusion protein of Amylase,alpha

亚型:IgG2b

克隆:360

描述:The complete human amylase gene contains five intact genes: two pancreatic genes (AMY2A and AMY2B) and three salivary genes (AMY1A,AMY1B,and AMY1C). AMY2A(Pancreatic alpha-amylase)is also named as PA and is a 56 kDa protein consisting of 496 amino acids in a single polypeptide chain,folded into three domains. It has a hydrolase activity acting on glycosyl bonds. AMY2A and AMY2B transcripts are present in the human pancreas but not in the parotid salivary gland,while AMY1 transcripts are present in the parotid gland but not in the pancreas.This antibody can recognize amylase alpha.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF