



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse ANXA2 Monoclonal Antibody

单克隆抗体

产品货号: MA01856

第三版

免疫原:Fusion protein of Annexin A2

亚型:IgG2b

克隆:99

描述:Annexin A2 (ANXA2),also named Annexin-2 or Lipocortin II,is a Ca²⁺ binding protein that is up-regulated in virally transformed cell lines and human tumors. This calcium-regulated membrane-binding protein,whose affinity for calcium is greatly enhanced by anionic phospholipids,may cross-link plasma membrane phospholipids with actin and the cytoskeleton and be involved with exocytosis. Heat-stress response may also involve Annexin-2.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF