

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ANP32B Monoclonal Antibody

单克隆抗体

产品货号: MA01857

第三版

免疫原:Fusion protein of ANP32B

亚型:IgG2b

克隆:44

描述:ANP32B is a member of the highly conserved acidic leucine-rich nuclear phosphoprotein 32 (ANP32) family. The ANP32 proteins have been implicated in a broad array of physiological processes,including cell differentiation,apoptotic cell death,and cell proliferation. ANP32B was identified as a novel direct substrate for caspase-3 and acted as a negative regulator for leukemic cell apoptosis. ANP32B controls the expression of the dendritic cell maturation factor CD83 by regulating the transport of its mRNA to the cytoplasm. It also modulates the activity of the transcription factor Kruppel-like factor 5 (KLF5).

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF