



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human APOD Monoclonal Antibody

单克隆抗体

产品货号: MA01860

第三版

免疫原:Fusion protein of APOD

亚型:IgG1

克隆:68

描述:Apolipoprotein D (ApoD) is a member of the lipocalin superfamily of ligand transporters, and has been implicated in the transport of small hydrophobic molecules. ApoD is also a component of plasma high-density lipoproteins (HDL). Alteration of ApoD expression has been linked to multiple neurological disorders, including Alzheimer's disease.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF