



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human APOC2 Monoclonal Antibody

单克隆抗体

产品货号: MA01862

第三版

免疫原:Fusion protein of Apolipoprotein H

亚型:IgG1

克隆:110

描述:Apolipoprotein H (ApoH),also known as beta2-glycoprotein I (B2-GPI),is a plasma glycoprotein,primarily synthesized in the liver. It has an important function in blood coagulation and clearance of apoptotic bodies from the circulation. ApoH is also an important actor of host innate immune response through its capacity to bind with high affinity to a large panel of pathogens or their proteins.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF