



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ARHGDIB Monoclonal Antibody

单克隆抗体

产品货号: MA01863

第三版

免疫原:Fusion protein of ARHGDIB

亚型:IgG2b

克隆:

描述:ARHGDIB,also known as RhoGDI2,is a conserved member of the RhoGDI family and plays an important role in cell migrations. It regulates the GDP/GTP exchange reaction of the Rho proteins by inhibiting the dissociation of GDP from them,and the subsequent binding of GTP to them. It is mainly in hematopoietic,endothelial,and epithelial cells. It has been linked to tumorigenesis and metastasis. RhoGDI2 expression is downregulated in several cancer types,such as bladder,lung and lymphoma,but is upregulated in prostate and gastric cancer.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF