



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ASS1 Monoclonal Antibody

单克隆抗体

产品货号: MA01864

第三版

免疫原:Fusion protein of ASS1

亚型:IgG1

克隆:

描述:ASS1 is a rate-limiting enzyme in the urea cycle. Citrullinemia type I is an autosomal recessive disorder that is caused by a deficiency of the urea cycle enzyme argininosuccinate synthetase (ASS1). Deficiency of ASS1 shows various clinical manifestations encompassing severely affected patients with fatal neonatal hyperammonemia as well as asymptomatic individuals with only a biochemical phenotype.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF