

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse ATG5 Monoclonal Antibody

单克隆抗体

产品货号: MA01866

第三版

免疫原:Fusion protein of ATG5

亚型:IgM

克隆:

描述:ATG5,also named as APG5L and ASP,belongs to the ATG5 family. It is required for autophagy. It plays an important role in the apoptotic process,possibly within the modified cytoskeleton. Its expression is a relatively late event in the apoptotic process,occurring downstream of caspase activity. Autophagy is a catabolic process for the autophagosomic-lysosomal degradation of bulk cytoplasmic contents. Formation of the autophagosome involves a ubiquitin-like conjugation system in which Atg12 is covalently bound to Atg5 and targeted to autophagosome vesicles. It mediates autophagosome-independent host protection.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF