

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ATP5C1 Monoclonal Antibody

单克隆抗体

产品货号: MA01867

第三版

免疫原:Fusion protein of ATP5C1

亚型:IgG2a

克隆:

描述:ATP5C1(ATP synthase subunit gamma,mitochondrial) is also named as ATP5C,ATP5CL1 and belongs to the ATPase gamma chain family. Some scientists reported the complete sequence of the gene for the human ATP synthase gamma subunit and described tissue-specific isoforms of the subunit generated by alternative splicing of exon 9. The liver (L) isoform differed from the heart (H) isoform by the addition of a single amino acid (asp273) at the C terminus. .

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF