



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human AZGP1 Monoclonal Antibody

单克隆抗体

产品货号: MA01868

第三版

免疫原:Fusion protein of AZGP1

亚型:IgG1

克隆:78

描述:Zinc-alpha-2-glycoprotein (AZGP1) is a 41-kDa soluble protein normally found in body fluids,functions as a lipid mobilizing factor . It is known to be expressed in the secretory epithelia of the liver,lung,breast,GI tract and sweat glands,sharing significant structural similarity with the class I major histocompatibility complex (MHC) antigens . AZGP1 is involved in carcinogenesis and differentiation. Altered expression of AZGP1 has been reported in breast cancer,prostate cancer and lung adenocarcinoma,hepatocellular carcinoma,pancreatic carcinoma and oral tumors .

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF