

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human B2M Monoclonal Antibody

单克隆抗体

产品货号: MA01869

第三版

免疫原:Synthetic Peptide

亚型:IgG2b

克隆:80

描述:Beta-2-microglobulin (B2M) is a component of MHC class I molecules, which are present on the surface of nearly all nucleated cells. It can be found in body fluids under physiologic conditions as a result of shedding from cell surfaces or intracellular release. B2M has various biological functions, including antigen presentation. Investigations reveal that increased synthesis and release of B2M are present in several malignant diseases.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF