

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat BAX Monoclonal Antibody

单克隆抗体

产品货号: MA01870

第三版

免疫原:Fusion protein of BAX

亚型:IgG2b

克隆:44

**描述:**BAX,also named as BCL2L4,is a pro-apoptotic member of the Bcl-2 protein family,which plays a pivotal role in controlling cell life and death. Bax largely localizes to the cytoplasm of healthy cells,but accumulates on the outer mitochondrial membrane upon apoptosis induction . BAX can commit a cell to apoptosis by translocation from the cytosol to the mitochondria and permeabilization of the outer mitochondrial membrane,which leads to the release of cytochrome c from mitochondria . The expression of BAX is upregulated by the tumor suppressor protein p53,and BAX has been shown to be involved in p53-mediated apoptosis .

**配方:**

**浓度:**None

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF