



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat BCL2 Monoclonal Antibody

单克隆抗体

产品货号: MA01872

第三版

免疫原:Fusion protein of BCL2

亚型:IgG2b

克隆:502

描述:BCL2 belongs to the Bcl-2 family. It suppresses apoptosis in a variety of cell systems including factor-dependent lymphohematopoietic and neural cells. BCL2 regulates cell death by controlling the mitochondrial membrane permeability. It appears to function in a feedback loop system with caspases. BCL2 inhibits caspase activity either by preventing the release of cytochrome c from the mitochondria and/or by binding to the apoptosis-activating factor (APAF-1).

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF