



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human BCL-x Monoclonal Antibody

单克隆抗体

产品货号: MA01873

第三版

免疫原:Fusion protein of Bcl-XL

亚型:IgG2b

克隆:

描述:BCL2L1 is a member of the BCL-2 protein family. BCL2L1 is expressed as two distinct isoforms, Bcl-xL and Bcl-xS, and is located at the outer mitochondrial membrane. The Bcl-XL isoform is a 233 amino acid protein, acting as an apoptotic inhibitor. Bcl-XL can form heterodimers with BAX, BAK or BCL2, and the heterodimerization with BAX does not seem to be required for anti-apoptotic activity. The Bcl-XS isoform is a shorter variant that is 178 amino acids in length and lacks a 63 amino acid region (amino acids 126-188), acting as an apoptotic activator.

配方:

浓度:None

保存: +4 °C 保存 (1-2 周). 长期保存在-20 °C 或者-70 °C. 避免反复冻融.

应用:ELISA, WB, ICC/IF