

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human BEST1 Monoclonal Antibody

单克隆抗体

产品货号: MA01875

第三版

免疫原:Fusion protein of BEST1

亚型:IgG1

克隆:58

描述:Bestrophin-1 (BEST1) is a 68 kDa transmembrane protein which belongs to the bestrophin family of anion channels. It is predominantly expressed in the basolateral membrane of the retinal pigment epithelium. Studies show that bestrophin 1 functions as a Ca(2+)-dependent Cl(-) channel and modulator of voltage-dependent L-type Ca(2+) channels. It may also contribute to the basolateral cell conductance in airway epithelial cells. This protein is encoded by the VMD2 gene. Mutations in VMD2 can lead to different types of retinal or macular degenerations, including Best vitelliform macular dystrophy (BMD), autosomal recessive bestrophinopathy (ARB), autosomal dominant vitreoretinopathy (ADVIRC) and adult-onset vitelliform macular dystrophy (AVMD).

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF