

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human BID Monoclonal Antibody

单克隆抗体

产品货号: MA01876

第三版

免疫原:Fusion protein of BID

亚型:IgG1

克隆:520

描述:BID,also named as p22 BID,can be cleaved into p11 BID,p13 BID and p15 BID. It is pro-apoptotic molecules. The major proteolytic product p15 BID allows the release of cytochrome c. BID-L,BID-EL and BID-ES induce ICE-like proteases and apoptosis. BID-S does not induce apoptosis. BID counters the protective effect of Bcl-2. .

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF