



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human BRD8 Monoclonal Antibody

单克隆抗体

产品货号: MA01877

第三版

免疫原:Fusion protein of BRD8

亚型:IgG1

克隆:407

**描述:**BRD8 is one of the bromodomain-containing proteins,which is an acetylated lysine- binding domain and thought to be involved in regulation of protein acetylation and/or HAT (histone acetyl transferase) activity . It is also reported to be an accessory subunit of NuA4 HAT complex (also known as TRRAP/TIP60 complex) through biochemical purification . In addition,it can bind to thyroid hormone receptor-  $\beta$  and/or retinoid X receptor,and can act as a coactivator of the nuclear hormone receptor-mediated transcription in reporter assays .

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF