

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Bovine BSA Monoclonal Antibody

单克隆抗体

产品货号: MA01879

第三版

免疫原:Synthetic Peptide

亚型:IgG1

克隆:22

描述:BSA, also named as ALB and Bovine Serum albumin, is the main protein of plasma. It has a good binding capacity for water, Ca²⁺, Na⁺, K⁺, fatty acids, hormones, bilirubin and drugs. Its main function is the regulation of the colloidal osmotic pressure of blood. This antibody recognizes albumin from bovine but not human, rat nor mouse.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF