

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human C1QC Monoclonal Antibody

单克隆抗体

产品货号: MA01880

第三版

免疫原:Fusion protein of C1QC

亚型:IgM

克隆:115

描述:The first component of complement,C1,is a calcium-dependent complex of the 3 subcomponents C1q,C1r,and C1s. C1q is composed of 18 polypeptide chains: six A-chains,six B-chains,and six C-chains. Each chain contains a collagen-like region located near the N terminus and a C-terminal globular region. Deficiency of C1q has been associated with lupus erythematosus and glomerulonephritis. This antibody is raised against C1qC which is the C-chain polypeptide of human complement subcomponent C1q.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF