

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human C4A (gamma chain) Monoclonal Antibody

单克隆抗体

产品货号: MA01881

第三版

免疫原:Fusion protein of C4

亚型:IgG1

克隆:30

描述: Complement C4, a component of the classical complement pathway, has an important role in innate immune function. C4 protein is expressed as a single chain precursor which is proteolytically cleaved into a trimer of alpha, beta, and gamma chains (molecular weights 95, 78, and 31 kDa) prior to secretion. The alpha chain is cleaved to release C4 anaphylatoxin, an antimicrobial peptide and a mediator of local inflammation. This antibody raised against 1543-1744aa of human C4A protein recognizes C4 gamma chain.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF