

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat C9orf72 Monoclonal Antibody

单克隆抗体

产品货号: MA01882

第三版

免疫原:Fusion protein of C9orf72

亚型:IgG2a

克隆:38

描述:C9ORF72 has a domain with polymorphic hexanucleotide repeat (GGGGCC). The C9ORF72-hexanucleotide repeat expansions have been recently identified as genetic markers in amyotrophic lateral sclerosis (ALS) and frontotemporal lobar degeneration (FTLD). FTLD-TDP has five subtypes: Sporadic FTLD, GRN mutation FTLD, TARDBP mutation FTLD, VCP mutation FTLD and C9ORF72 mutation FTLD. The C9ORF72 repeat expansions may indicate a worse prognosis in ALS. Human C9ORF72 has some isoforms with MW 54-60 kDa and 25-30 kDa. Mouse C9orf72 has some isoforms with MW 50-60 kDa and 35 kDa.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF