



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CAPG Monoclonal Antibody

单克隆抗体

产品货号: MA01883

第三版

免疫原:Fusion protein of CAPG

亚型:IgG1

克隆:40

**描述:**CapG (gelsolin-like actin-capping protein) is a member of the gelsolin superfamily proteins, which regulates actin filament length by capping or severing filaments. Studies have shown a role for CapG in regulating cell motility. The expression of CapG is upregulated under hypoxia and in turn increases the motility of cells. Recently it has been reported that the overexpression of CapG is associated with the progression of lung adenocarcinoma as well as lymphatic invasion status of cholangiocarcinoma, which may be used for prognosis of these tumors. (21908955, 22155129)

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF