

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse CASP8 Monoclonal Antibody

单克隆抗体

产品货号: MA01885

第三版

免疫原:Fusion protein of Caspase 8

亚型:IgG2b

克隆:98

描述:Caspase 8,also named as MCH5,MACH,FLICE,and CAP4,belongs to the peptidase C14A family. It may participate in the GZMB apoptotic pathways and functions as an upstream regulator in a-bisabolol-induced apoptosis. Caspase 8 catalyzes an essential intermediate step in the ubiquitination and proteasome-mediated degradation of IRF3. It may control diabetic embryopathy-associated apoptosis via Caspase 8 is expressed as nine isoforms by alternative splicing with the molecular mass from 26 kDa to 62 kDa. Caspase 8 Monoclonal Antibody can recognize pro- and cleaved-caspase 8.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF