



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CAV1 Monoclonal Antibody

单克隆抗体

产品货号: MA01887

第三版

免疫原:Fusion protein of Caveolin-1

亚型:IgG1

克隆:60

**描述:**Caveolin-1 (CAV1),a multifunctional protein,is the main constituent molecule of caveolae and represents a scaffolding molecule for several signaling molecules including epidermal growth factor receptor . Several studies have implicated that a reduced expression of CAV1 was found in cancers including head and neck carcinoma . However,other studies recognize CAV1 as a tumor promoter because CAV1 is overexpressed in various kinds of cancers,especially in oral cancer . Recent study also show that CAV1 is involved in astric Cancer .

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF