

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD13 Monoclonal Antibody

单克隆抗体

产品货号: MA01890

第三版

免疫原:Fusion protein of CD13

亚型:IgG2a

克隆:339

描述:CD13/aminopeptidase N is a 150 kDa cell surface glycoprotein,which was originally used as a marker for subpopulations of hematopoietic cells. CD13 is widely distributed in a variety of mammalian cells,such as monocyte/macrophages, fibroblasts,neutrophils,endothelial cells,and epithelial cells. CD13 is also a pan-myeloid marker,present on mature granulocytes and monocytes.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF