



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD138 Monoclonal Antibody

单克隆抗体

产品货号: MA01892

第三版

免疫原:Fusion protein of CD138/Syndecan-1

亚型:IgG2a

克隆:666

描述:Syndecan-1 (SDC1),also named CD138, is an integral membrane protein. It participates in cell proliferation,cell migration and cell-matrix interactions via its receptor for extracellular matrix proteins. It is a heparan sulfate proteoglycan expressed on the surface of,and actively shed by,myeloma cells. Altered syndecan-1 expression has been detected in several different tumor types. SDC1 was regarded as a useful marker for labeling normal and neoplastic plasma cells and plasmacytoid lymphomas.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF