

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD1A Monoclonal Antibody

单克隆抗体

产品货号: MA01895

第三版

免疫原:Fusion protein of CD1A

亚型:IgG2b

克隆:479

描述:CD1A is an antigen-presenting protein that binds self and non-self lipid and glycolipid antigens and presents them to T-cell receptors on natural killer T-cells. During protein synthesis and maturation,CD1 family members bind endogenous lipids that are replaced by lipid or glycolipid antigens when the proteins are internalized and pass through endosomes,before trafficking back to the cell surface. CD1A is a 50-kDa glycoprotein. It has a heterodimer form with molecular weight of 65-70 kDa.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF