



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD20 Monoclonal Antibody

单克隆抗体

产品货号: MA01896

第三版

免疫原:Fusion protein of CD20

亚型:IgG1

克隆:103

描述:CD20 is a 33-37 kDa transmembrane phosphoprotein. CD20 is a B-lymphocyte surface molecule that is widely expressed during B-cell ontogeny, from early pre-B-cell developmental stages until final differentiation into plasma cells. CD20 functions as calcium-permeable cation channel. It is involved in the regulation of B-cell activation and proliferation. CD20 serves as a useful target for antibody-mediated therapeutic depletion of B-cells.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF