



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse CD23 Monoclonal Antibody

单克隆抗体

产品货号: MA01899

第三版

免疫原:Fusion protein of CD23

亚型:IgG2a

克隆:168

描述:CD23,also known as low affinity immunoglobulin epsilon Fc receptor,is a transmembrane glycoprotein present on a subpopulation of B lymphocytes in germinal centres,EBV-transformed B-lymphoblastoid cell lines,follicular dendritic cells,and a subpopulation of peripheral blood cells. CD23 has essential roles in the regulation of IgE production and in the differentiation of B-cells.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF