
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD43 Monoclonal Antibody

单克隆抗体

产品货号: MA01905

第三版

免疫原:Fusion protein of CD43

亚型:IgG2b

克隆:46

描述:CD43,also named leukosialin,sialophorin,galactoglycoprotein,leukocyte sialoglycoprotein,SPN,and is a mucin-like type I transmembrane protein. In humans,it is expressed in haematopoietic cells,including T lymphocytes,monocytes,granulocytes,natural killer cells,platelets,and haematopoietic stem cells,with the exception of mature erythrocytes and B cell subpopulations. The unglycosylated form of CD43 migrates on SDS-PAGE with a molecular weight of 54 kDa ,whereas CD43 from different haematopoietic cell lines displays increased molecular weights since it is glycosylated with O-linked chains that differ in core structure and sialylation.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF