

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD44 Monoclonal Antibody

单克隆抗体

产品货号: MA01906

第三版

免疫原:Fusion protein of CD44

亚型:IgG2a

克隆:23

描述:CD44 is a type I transmembrane glycoprotein that mediates cell-cell and cell-matrix interactions through its affinity for hyaluronic acid (HA) and possibly also through its affinity for other ligands. Adhesion with HA plays an important role in cell migration,tumor growth and progression. CD44 is also involved in lymphocyte activation,recirculation and homing,and in hematopoiesis. This protein exists in multiple isoforms (with molecular weight ranging from 16 kDa to over 200 kDa) generated by alternative RNA splicing and undergoes extensive post-translational modifications.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF