

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat CD9 Monoclonal Antibody

单克隆抗体

产品货号: MA01911

第三版

免疫原:Fusion protein of CD9

亚型:IgG1

克隆:129

**描述:**The cell-surface molecule CD9,a member of the transmembrane-4 superfamily,interacts with the integrin family and other membrane proteins,and is postulated to participate in cell migration and adhesion. Expression of CD9 enhances membrane fusion between muscle cells and promotes viral infection in some cells . It is often used as a mesenchymal stem cell marker . CD9 is also known as the p24 antigen besides MIC3,TSPAN29 because it is a protein of molecular weight 24 kD. The CD9 antigen appears to be a 227-amino acid molecule with 4 hydrophobic domains and 1 N-glycosylation site.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF