

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CDX2 Monoclonal Antibody

单克隆抗体

产品货号: MA01914

第三版

免疫原:Fusion protein of CDX2

亚型:IgG1

克隆:913

**描述:**CDX2,also named as Homeobox protein CDX-2,is a 313 amino acid protein,which contains one homeobox DNA-binding domain and belongs to the Caudal homeobox family. CDX2 localizes in the nucleus and is involved in the transcriptional regulation of multiple genes expression in the intestinal epithelium. The relative expression of CDX1 to CDX2 protein may be important in the anterior to posterior patterning of the intestinal epithelium and in defining patterns of proliferation and differentiation along the crypt-villus axis. Both Cdx1 and Cdx2 genes must be expressed to reduce tumorigenic potential,to increase sensitivity to apoptosis,and to reduce cell migration,suggesting that the two genes control the normal phenotype by independent pathways.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF