



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CP Monoclonal Antibody

单克隆抗体

产品货号: MA01915

第三版

免疫原:Fusion protein of Ceruloplasmin

亚型:IgG1

克隆:11B

描述:Ceruloplasmin (CP) is a serum glycoprotein synthesized by hepatocytes and functions as the major transport protein of copper in blood where it acts as ferroxidase and is required for iron transport of transferrin. The intact ceruloplasmin migrated around 130-140 kDa and some proteolytic cleaved fragments could also be observed.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF