



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CHEK1 Monoclonal Antibody

单克隆抗体

产品货号: MA01916

第三版

免疫原:Fusion protein of Chk1

亚型:IgG2b

克隆:32

**描述:**CHEK1(Checkpoint kinase-1) is also named as CHK1 and belongs to the protein kinase superfamily. It is implicated in a circuit in which it activates checkpoints,DNA repair and proliferating cell nuclear antigen and FANCD2 monoubiquitinylation. CHEK1 protects vertebrate cells against spontaneous chromosome missegregation and is required to sustain anaphase delay when spindle function is disrupted by taxol. It has 3 isoforms produced by alternative splicing with the molecular mass of 54 kDa,44 kDa and 50 kDa.

**配方:**

**浓度:**None

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF