



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CHRNA7 Monoclonal Antibody

单克隆抗体

产品货号: MA01918

第三版

免疫原:Fusion protein of CHRNA7

亚型:IgG2b

克隆:34

描述:Nicotinic acetylcholine receptors (nAChRs) are cholinergic receptors that form ligand-gated ion channels that mediate fast signal transmission at synapses. CHRNA7 (neuronal acetylcholine receptor subunit alpha-7) forms a homo-oligomeric channel, displays marked permeability to calcium ions and is a major component of brain nicotinic receptors that are blocked by, and highly sensitive to, alpha-bungarotoxin.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF