

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CLEC11A Monoclonal Antibody

单克隆抗体

产品货号: MA01922

第三版

免疫原:Fusion protein of CLEC11A

亚型:IgG2b

克隆:239

描述:CLEC11A,also named as CLECSF3,LSLCL,SCGF and p47,stimulates the proliferation and differentiation of hematopoietic precursor cells from various lineages,including erythrocytes,lymphocytes,granulocytes and macrophages. CLEC11A acts synergistically with other cytokines,including IL-3,GCSF,GMCSF and FLT3 ligand. It suppresses SCF-stimulated erythrocyte proliferation.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF