



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CLTCL1 Monoclonal Antibody

单克隆抗体

产品货号: MA01924

第三版

免疫原:Synthetic Peptide

亚型:IgG2a

克隆:55

描述:CLTCL1,also named as CLH22,CLTCL and CLTD,belongs to the clathrin heavy chain family. Clathrin is the major protein of the polyhedral coat of coated pits and vesicles. Two different adapter protein complexes link the clathrin lattice either to the plasma membrane or to the trans-Golgi network. The antibody is specific to CLTCL1. It has no cross reaction to CLTC.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF