



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CNN2 Monoclonal Antibody

单克隆抗体

产品货号: MA01925

第三版

免疫原:Fusion protein of CNN2

亚型:IgG1

克隆:33

描述:Calponin is a family of actin filament-associated proteins which regulate smooth muscle cell contraction. Three isoforms of calponin exist: calponin h1 (CNN1), calponin h2 (CNN2) and calponin 3 (CNN3). Calponin 1 and calponin 2 are predominately expressed in smooth muscle cells and cardiac muscle cells, respectively. Calponin 3 is highly expressed in many tissues including articular cartilage and brain.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF