



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CNPY2 Monoclonal Antibody

单克隆抗体

产品货号: MA01926

第三版

免疫原:Fusion protein of CNPY2/ MSAP

亚型:IgG2b

克隆:448

描述:CNPY2,also named as MSAP, TMEM4 and ZSIG9,is a positive regulator of neurite outgrowth. It interacts with the ezrin-radixin-moesin (ERM)-like myosin regulatory light chain-interacting protein (MIR),and the two proteins are co-localized in cell lines and in primary neurons. CNPY2 can stabilize MRLC and thus bring about an increase in neurite outgrowth. This antibody is specific to CNPY2

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF