



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CFB Monoclonal Antibody

单克隆抗体

产品货号: MA01928

第三版

免疫原:Fusion protein of Complement factor B

亚型:IgM

克隆:12

描述:Complement factor B (CFB) is a component of the alternative pathway of complement activation. Complement factor B is cleaved into Ba and Bb by factor D in the presence of C3b. The active subunit Bb is a serine protease which associates with C3b to form the alternative pathway C3 convertase. Bb is involved in the proliferation of preactivated B lymphocytes, while Ba inhibits their proliferation.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF